



PentaBase

PRODUCT CATALOGUE

Customised
Oligonucleotides
and PCR Assays

2024

Contact

PentaBase A/S
Petersmindevej 1A
DK-5000 Odense C

+45 36 96 94 96

Info@pentabase.com

www.pentabase.com

PentaBase

Inventors of INA[®] - Experts in real-time PCR

PentaBase is a Danish science-driven company established in 2007 by CEO Ulf Bech Christensen. Our products and services are based on unique DNA technologies providing innovative and superior PCR workflows. Our mission is to provide PCR-based assays that enable personalised treatment using sensitive, specific and robust monitoring of genetic biomarkers.

Our approach is to continuously innovate and improve our products and services to provide novel or improved PCR-based solutions for the benefit of patients.

We see every patient as a unique individual who deserves to be treated individually. It is pretty simple: We wish to participate in creating a healthier society through personalised detection, prevention, and treatment of diseases, based on precise and timely genetic diagnostics.

We strive to use our expertise in PCR to offer our customers the best service possible and we aim to combine fast delivery of superior quality products with optimal support before, during, and after the purchase.

Specifications subject to change without notice.
Copyright PentaBase A/S.

Version: 2.5
Last revised: November 2023







	INA®	04
	Probes	05
	EasyBeacon™	05
	HydroEasy®	06
	Fluorophores & Quenchers	07
	BaseBlockers™ and SuPrimers	08
	BaseBlockers™	08
	SuPrimers™	09
	Custom PCR Assays	12
	How to order	13
	Prices	14
	Oligonucleotides	14
	Assay design	16
	Custom individual services	17





Illustration created with BioRender.com

Intercalating Nucleic Acid - INA®

We supply customized oligonucleotides and PCR assays that are based on our unique and proprietary Intercalating Nucleic Acid (INA®) technology*. Intercalating Nucleic Acids are made by adding Intercalating Pseudo Nucleotide (IPN) modifications to standard DNA oligonucleotides during synthesis. IPNs are hydrophobic nucleoside analogues that are added to the growing oligonucleotide chain without replacing existing nucleotides.


INA® is the only DNA platform technology that works by increasing the π stacking effect of nucleotide base pairs. The result is a much higher affinity and specificity of INA® oligonucleotides towards the target nucleotide sequence compared to standard DNA oligonucleotides. This allows for enhanced discrimination against mismatches. In addition, temperature independent interactions between IPN modifications of non-target bound INA®-based fluorescent probes result in low background fluorescence at all temperatures.

Our EasyBeacon™ and HydrolEasy® probes form the backbone of our fluorescence-based assays providing nuclease-resistance, high target specificity and high signal-to-noise ratio. Our SuPrimers™ reduce primer-dimer formation and reduce off-target binding, and our BaseBlockers™ allow for extremely selective amplification of mutated DNA in a wild type DNA background.

At PentaBase A/S we have developed a wide range of INA®-based PCR applications for our customers as well as for our in vitro diagnostic assays for molecular oncology and infectious disease.


* 1. Christensen UB, Pedersen EB. 2002. <https://www.doi.org/10.1093/nar/gkf624>
2. Christensen UB, et al. 2004. <https://www.doi.org/10.1081/ncn-120027829>
3. Nielsen CB, et al. 2004. Doi: <https://www.doi.org/10.1021/bc0341932>

EasyBeacon™ probes



The perfect probe for end-point melt analyses. The INA® technology used in EasyBeacon™ probes creates a secondary structure in the unbound state that efficiently quenches the signal resulting in a very low background fluorescence without adding a stem sequence. The IPNs in the probes furthermore prevent hydrolysis, making the probes nuclease resistant and therefore intact after amplification.

The high specificity and low background makes it possible to easily differentiate between single nucleotide polymorphisms (SNPs) as well as multiplexing in end-point melting curve analyses.




The example shows an EasyBeacon™ assay analysing a human genomic DNA target region with a single nucleotide polymorphism (SNP). A difference in melting temperature of 4°C between the wild type target and the SNP target makes it easy to distinguish between homozygous wild type (blue line), homozygous mutant (green line) and heterozygous (red line).


EasyBeacon™ probe advantages

- Ideal for end-point melt analyses.
- Nuclease resistant – intact after amplification.
- Unbound probe is quenched at all temperatures - no need to add a stem
- Ideal for multiplexed genotyping and High Resolution Melt (HRM) applications

HydroEasy[®] probes



The HydroEasy[®] probe is an improved alternative to TaqMan hydrolysis probes. The INA[®] technology used in HydroEasy[®] probes creates a secondary structure in the unbound state that efficiently quenches the signal resulting in a very low background fluorescence and consequently higher signal-to-noise ratio compared to standard hydrolysis probes.



The example shows the raw fluorescence and normalized fluorescence data of two HydroEasy[®] probe assays compared to standard hydrolysis probe assay. The lower background fluorescence of the HydroEasy[®] probes results in significantly higher signal-to-noise (S/N) ratios.

HydroEasy[®] probe advantages

- Unbound probe is quenched at all temperatures
- Enhanced signal-to-noise ratio
- Improved target specificity

Fluorophores & Quenchers

	Fluorophore	Excitation	Emission	Quenchers				Quencher	Range
	ATTO 390	390 nm	476 nm	BHQ-1 λ (max)=534nm Range=480-580nm	Dabcyl λ (max)=453nm Range=380-530nm	Tamra λ (max)=544nm Range=470-560nm	Penta GQ/YQ λ (max)=531nm Range=480-580nm	BHQ™-1	480-580 nm
	ATTO 425	439 nm	485 nm					BHQ™-2	550-650 nm
	ATTO 465	453 nm	509 nm					BBQ-650™	550-750 nm
	Penta Green	494 nm	518 nm					Dabcyl	380-530 nm
	6-FAM	494 nm	518 nm					Tamra	470-560 nm
	ATTO 495	498 nm	526 nm					Penta GQ/YQ	480-580 nm
	ATTO 490LS	498 nm	658 nm						
	ATTO 488	500 nm	520 nm						
	Tide Fluor™ 2 (TF2)	503 nm	525 nm						
	6-JOE	520 nm	548 nm						
	6-TET	521 nm	536 nm						
	CAL Fluor® Gold 540	522 nm	543 nm						
	ATTO 532	532 nm	552 nm						
	6-HEX	535 nm	556 nm						
	CAL Fluor® Orange 560	537 nm	558 nm						
	Penta Yellow	538 nm	551 nm						
	Quasar® 570	547 nm	570 nm						
	CY-3	550 nm	570 nm						
	Penta Dark Yellow	550 nm	570 nm						
	ATTO 550	554 nm	576 nm						
	Tamra	555 nm	580 nm						
	ATTO 565	564 nm	590 nm						
	Penta Orange	577 nm	601 nm						
	ROX	578 nm	602 nm						
	Penta Light Orange	581 nm	594 nm						
	CAL Fluor® Red 610	590 nm	610 nm						
	California Red™	591 nm	608 nm						
	ATTO 590	593 nm	622 nm						
	ATTO 594	603 nm	626 nm						
	Quasar® 670	644 nm	670 nm						
	ATTO 647N	646 nm	664 nm						
	CY-5	650 nm	670 nm						
Penta Red	650 nm	670 nm							
Penta Purple	675 nm	694 nm							
Quasar® 705	690 nm	705 nm							
ATTO 700	700 nm	716 nm							
ATTO 725	728 nm	751 nm							

Quencher	λ (max)
BHQ™-1	534 nm
BHQ™-2	579 nm
BBQ-650™	650 nm
Dabcyl	453 nm
Tamra	544 nm
Penta GQ/YQ	531 nm

BaseBlocker™





Illustration created with BioRender.com

BaseBlocker™s are short oligonucleotides that are highly modified with several IPNs. BaseBlocker™s have high target affinity and specificity, allowing sequence specific blocking of for example wild type templates*. In addition, BaseBlockers™ are modified to ensure that they do not work as potential primers.

BaseBlocker™s enable detection of down to a single DNA molecule with a single SNP in a background of wild-type DNA[^], lowering the risk of false positive signals without influencing the rate of false negatives. The example shows how adding a BaseBlocker™ to an allele-specific qPCR Assay blocks amplification of the wild type allele (red and blue lines) without affecting amplification of the mutant allele (green and black lines).



BaseBlocker™ advantages

- INA® modified BaseBlocker™ for highly specific sequence discrimination
- Extremely sensitive mutant detection in wild-type background reducing the risk of false positives

* 1. Riva A, et al. 2017. <https://www.doi.org/10.1371/journal.pone.0178027>


2. Sørensen AL, et al. 2023. <https://www.doi.org/10.1371/journal.pone.0281558>

[^] Jensen SG, et al. 2021. <https://www.doi.org/10.1371/journal.pone.0253687>


SuPrimer™

The INA® modified SuPrimer™ works as a normal primer in the PCR reaction, just better. It gives you the capability of avoiding primer-dimer formation and at the same time increasing affinity, specificity, sensitivity and signal strength.

Amplification using regular primers often results in unspecific products due to primer-dimer interactions and/or other off target events. The increased affinity and specificity of the SuPrimer™s greatly reduces the risk of these primer-dimer interactions and off-target events, resulting in amplification of only the target sequence.




The example above shows an end-point melting curve following amplification of the RNase P gene used as internal control in SARS-CoV-2 testing with the primers and SuPrimer™s both made on the basis of the WHO recommendations. The assay based on SuPrimer™s show a only single PCR product compared to several products when using standard DNA primers and also higher signal strength.



Including the SuPrimer™s in qPCR assay design can increase specificity, sensitivity and signal strength. The example above shows the difference between the amplification curves of two normal primer sets compared to two corresponding primer sets comprising a modified reverse SuPrimer™. Primer1 and SuPrimer™1 amplifies an identical sequence and the same applies to Primer2 and SuPrimer™2. The use of SuPrimer™s results in higher sensitivity and higher fluorescent signals than the use of regular primers.

SuPrimer™ advantages

- Reduces primer-dimer formation
- Lower off-target binding and higher target specificity and affinity
- Higher signal strength




EpiDirect®

The EpiDirect® technology utilizes a type of INA® known as an EpiPrimer™ that allows for direct evaluation of DNA methylation by qPCR without prior chemical conversion (reference til Nature Genetics artiklen).

The EpiPrimer™ consists of 3 parts: a) an IPN-containing anchor sequence that binds to the methylated target region, b) a starter sequence that primes DNA replication, and c) a loop sequence that is part of the subsequent methyl-unspecific priming. Thus, the anchor and starter sequence prime the initial replication of the methylated DNA while the loop and starter sequence combined work as a primer in the remaining PCR cycles.

The high specificity of the IPN-containing anchor sequence towards the methylated target allows for highly selective amplification of the methylated sample DNA.



EpiDirect® advantages

- Direct evaluation of DNA methylation using customised INA®-based EpiPrimer™ design
- Avoids artefacts related to chemical conversion of DNA prior to real-time PCR



Custom PCR assays

As leading PCR experts, we have vast experience in developing and producing highly sensitive and specific assays for our customers. Our knowhow and unique INA™ technology allows us develop novel high performance PCR assays or significantly improve the performance of client assays previously based on standard nucleotide chemistries.

For our customized assays we have specific minimum quality standards for delivery and guaranteed performance:

- Minimum 10 nmol delivered
- Validation will be performed on PentaBase Mastermix and PentaBase BasePurifier™ / BaseTyper™ instruments
- SNP genotyping assay (melt-based) will find 100 or less copies of synthetic SNP variant target in 5 ng synthetic wild-type background with a specificity determined by a temperature difference between variant and wild-type target of minimum 4 degrees
- Pathogen identification assay with a limit of detection (LOD) of 100 or less copies of synthetic template hidden in a background of 5 ng human genomic DNA
- Somatic mutation assays with a LOD of 1% or better measured using synthetic mutated templates hidden in a background of human genomic DNA
- Custom EpiDirect® Anchor Design specificity as determined by a difference in melting temperature between methylated and unmethylated target of at least 5°C
- Custom EpiDirect Assay Design limit of detection (LoD) of 5% methylation or better

How to order

Customised oligonucleotides and PCR assays can be ordered by sending an email to **order@pentabase.com**. When ordering custom oligonucleotides, please include the following information, and do not worry about including the IPNs in your design, our specialists will add IPNs once a query has been sent:

- Synthesis scale
- Sequence (5' to 3') + given sequence name
- 5' and 3' modifications
- Purification
- Technology

Synthesis scale	Technology	Purification
10 nmol delivered	EasyBeacon™ probe	RPLC
50 nmol delivered	HydroEasy® probe	HPLC
100+ nmol delivered (discounts apply)	SuPrimer™	2x HPLC
	BaseBlocker™	
	EpiPrimer™	
	Conventional primer	
	Conventional probe	

Our custom oligonucleotides can be ordered in fixed amounts from 10 nmol delivered and supplied as either lyophilized or in aliquoted solutions as specified by the customer. We use several control steps to ensure oligonucleotide quality including measuring trityl data to calculate synthesis efficiency, evaluating HPLC spectra to ensure pure oligonucleotide output, and measuring optical densities for accurate concentration determination. In addition, we perform a target melting curve analysis of our probes to evaluate target binding properties.

Purification can be done through RPLC (cartridge purification) or RP-HPLC, and every order is sent with a synthesis report. A copy of the melt analysis can be included upon request.

Prices

Below is the price table for our custom oligonucleotides.

DNA Primer™ 14-39 mer		
Ref. No.	Product name	List price
430	DNA Primer; 10 nmol delivered; RPLC	50,00 €
431	DNA Primer; 10 nmol delivered; HPLC	60,00 €
432	DNA Primer; 50 nmol delivered; RPLC	60,00 €
433	DNA Primer; 50 nmol delivered; HPLC	70,00 €
DNA Primer™ 40-100 mer		
Ref. No.	Product name	List price
690	DNA Primer; 10 nmol delivered; HPLC	120,00 €
691	DNA Primer; 50 nmol delivered; HPLC	140,00 €
SuPrimer™ 14-35 mer		
Ref. No.	Product name	List price
434	SuPrimer™; 10 nmol delivered; HPLC	90,00 €
435	SuPrimer™; 50 nmol delivered; HPLC	100,00 €
BaseBlocker™ 14-35 mer		
Ref. No.	Product name	List price
436	Design of BaseBlocker™	500,00 €
437	BaseBlocker™; 10 nmol delivered; HPLC	140,00 €
438	BaseBlocker™; 50 nmol delivered; HPLC	280,00 €
EasyBeacon™ / HydrolEasy® Probe 14-35 mer; 10 nmol Delivered		
Ref. No.	Product name	List price
439/459	FAM/PentaGreen-BHQ1/Green Quencher/DPQ1	160,00 €
440/460	HEX/PentaYellow-BHQ1/Yellow Quencher/DPQ1	180,00 €
441/461	CalFluor Red/PentaOrange-BHQ2/ODPQ2	230,00 €
442/462	Atto 425/465-BHQ1/DPQ1	230,00 €
443/463	Atto 565-BHQ2/DPQ2	230,00 €
444/464	Quasar 670/Cy5/PentaRed-BHQ2/DPQ2	300,00 €
445/465	Atto 390/495-BHQ1/DPQ1	300,00 €
446/466	Quasar 705/Cy5.5/PentaPurple-BHQ2/DPQ2	350,00 €
447/467	Atto 488/490/532/550-BHQ1/DPQ1	460,00 €
448/468	Atto 594/647N/700/725-BHQ2/BBQ-650/DPQ2	460,00 €

EasyBeacon™ / HydrolEasy® Probe 14-35 mer; 50 nmol Delivered		
Ref. No.	Product name	List price
449/469	FAM/PentaGreen-BHQ1/Green Quencher/DPQ1	240,00 €
450/470	HEX/PentaYellow-BHQ1/Yellow Quencher/DPQ1	300,00 €
451/471	CalFluor Red/PentaOrange-BHQ2/ODPQ2	480,00 €
452/472	Atto 425/465-BHQ1/DPQ1	490,00 €
453/473	Atto 565-BHQ2/DPQ2	650,00 €
454/474	Quasar 670/Cy5/PentaRed-BHQ2/DPQ2	650,00 €
455/475	Atto 390/495-BHQ1/DPQ1	660,00 €
456/476	Quasar 705/Cy5.5/PentaPurple-BHQ2/DPQ2	700,00 €
457/477	Atto 488/490/532/550-BHQ1/DPQ1	920,00 €
458/478	Atto 594/647N/700/725-BHQ2/BBQ-650/DPQ2	920,00 €
DNA Probe 14-35 mer; 10 nmol Delivered		
Ref. No.	Product name	List price
479	FAM/PentaGreen-BHQ1/Green Quencher/DPQ1	130,00 €
480	HEX/PentaYellow-BHQ1/Yellow Quencher/DPQ1	150,00 €
481	CalFluor Red/PentaOrange-BHQ2/ODPQ2	200,00 €
482	Atto 425/465-BHQ1/DPQ1	200,00 €
483	Atto 565-BHQ2/DPQ2	200,00 €
484	Quasar 670/Cy5/PentaRed-BHQ2/DPQ2	270,00 €
485	Atto 390/495-BHQ1/DPQ1	270,00 €
486	Quasar 705/Cy5.5/PentaPurple-BHQ2/DPQ2	320,00 €
487	Atto 488/490/532/550-BHQ1/DPQ1	430,00 €
488	Atto 594/647N/700/725-BHQ2/BBQ-650/DPQ2	430,00 €
DNA Probe 14-35 mer; 50 nmol Delivered		
Ref. No.	Product name	List price
489	FAM/PentaGreen-BHQ1/Green Quencher/DPQ1	180,00 €
490	HEX/PentaYellow-BHQ1/Yellow Quencher/DPQ1	240,00 €
491	CalFluor Red/PentaOrange-BHQ2/ODPQ2	420,00 €
492	Atto 425/465-BHQ1/DPQ1	430,00 €
493	Atto 565-BHQ2/DPQ2	430,00 €
494	Quasar 670/Cy5/PentaRed-BHQ2/DPQ2	590,00 €
495	Atto 390/495-BHQ1/DPQ1	600,00 €
496	Quasar 705/Cy5.5/PentaPurple-BHQ2/DPQ2	640,00 €
497	Atto 488/490/532/550-BHQ1/DPQ1	860,00 €
498	Atto 594/647N/700/725-BHQ2/BBQ-650/DPQ2	860,00 €

Discounts - Oligonucleotides*	
nmol delivered	Discount rate
100 nmol	10%
250 nmol	20%
500 nmol	25%
1000 nmol	30%
2000 nmol	35%

*Discounts are applied to 50 nmol delivered based on order volume

Below is the price table for our New Assay Services.

NEW SNP ASSAY (one base change)		
Ref. No.	Product name	List price
620	New SNP genotyping simplex assay design	2.460,00 €
621	New SNP genotyping duplex assay design	5.730,00 €
622	New SNP genotyping triplex assay design	8.990,00 €
623	New SNP genotyping quadruplex assay design	11.860,00 €
625	New SNP genotyping simplex assay verification	3.620,00 €
626	New SNP genotyping duplex assay verification	8.660,00 €
627	New SNP genotyping triplex assay verification	13.690,00 €
628	New SNP genotyping quadruplex assay verification	18.120,00 €
	New SNP genotyping assay custom services (specific performance requirements, performance evaluation, etc.)	<i>Inquire</i>
NEW PATHOGEN IDENTIFICATION ASSAY		
Ref. No.	Product name	List price
635	New simplex pathogen identification assay design	2.460,00 €
636	New duplex pathogen identification assay design	5.730,00 €
637	New triplex pathogen identification assay design	8.990,00 €
638	New quadruplex pathogen identification assay design	11.860,00 €
640	New simplex pathogen identification assay verification	3.620,00 €
641	New duplex pathogen identification assay verification	8.660,00 €
642	New triplex pathogen identification assay verification	13.690,00 €
643	New quadruplex pathogen identification assay verification	18.120,00 €
	New pathogen identification assay custom services (specific performance requirements, performance evaluation, etc.)	<i>Inquire</i>
NEW SOMATIC MUTATION ASSAY		
Ref. No.	Product name	List price
575	New simplex somatic mutation assay design	3.020,00 €
576	New duplex somatic mutation assay design	6.850,00 €
577	New triplex somatic mutation assay design	10.880,00 €
578	New quadruplex somatic mutation assay design	14.510,00 €
580	New simplex somatic mutation assay verification	3.620,00 €
581	New duplex somatic mutation assay verification	8.660,00 €
582	New triplex somatic mutation assay verification	13.690,00 €
583	New quadruplex somatic mutation assay verification	18.120,00 €
	New somatic mutation assay custom services (specific performance requirements*, performance evaluation, etc.)	<i>Inquire</i>

CUSTOM EpiDirect® DESIGN		
Ref. No.	Product name	List price
660	Custom EpiDirect® Anchor Design and Validation	3,630.00 €
661	Custom EpiDirect® Assay Design and Validation	29,330.00 €

Below is the price table for our Custom Individual Services.

CUSTOM INDIVIDUAL SERVICES		
Ref. No.	Product name	List price
	Custom assay analytical performance evaluation IVDR Stage 4 - Accuracy ¹ , Analytical Sensitivity ² and Analytical Specificity ³	<i>Inquire</i>
	Custom EpiDirect® Assay Optimisation ⁴	<i>Inquire</i>
	<p>¹Accuracy including <i>Trueness</i> and <i>Precision</i></p> <p><i>Trueness</i> will be determined by calculating limits of agreement (LOA) using relevant available certified reference materials and at least 20 measurement points.</p> <p><i>Precision</i> will include repeatability (n≥ 20) and within laboratory reproducibility for up to 3 factors (lots, time, instruments, etc.)</p> <p>²Analytical sensitivity will include limit of blank (LoB), limit of detection (Lod) and limit of quantitation (LoQ) when relevant, of at least 20 measurement points.</p> <p>³Analytical specificity will include relevant in silico analyses and wet lab testing of up to 5 interfering targets or substances.</p> <p>⁴Oligonucleotide titration, PCR program optimization, etc.</p> <p>All analyses require availability of clinically relevant sample material (FFPE, plasma, etc.) and related costs are not included.</p> <p>Data will be summarised in an analytical performance evaluation report</p>	





Scan this QR-code to download our instrument and extraction kit catalogue!



Scan this QR-code to download our PCR assay catalogue!



Scan this QR-code to download the latest digital copy of this catalogue!

PentaBase

PentaBase A/S
Petersmindevej 1A
DK-5000 Odense C
Phone +45 36 96 94 96

For further information, visit
www.pentabase.com